



# **Aesthetic, Plastic, Reconstructive Surgery and Dermatology**

Specialized Operating Tables and Chairs

# Welcome to BRUMABA

## Reliable quality for your health

**BRUMABA is a family-run company that specializes in the development and production of ergonomic operating tables, chairs and stools:**

**By developing high-quality, innovative products, we aim to improve the health and well-being of healthcare professionals and patients around the world.**

The principle of our work is sustainability: BRUMABA products have been manufactured in our Bavarian homeland since 1980. The use of high-quality materials and reliable technology provide a proven long service life of the products.

Convince yourself of our noticeable advantage in ergonomics, comfort, reliability and efficiency.

Reliability that counts. Stability that lasts.  
BRUMABA - your partner in quality.

**MADE IN GERMANY**



**SINCE 1980**



Operating chair  
**BALANCE ADVANCE**

Patient recovery chair  
**CADDY**



Operating table  
**VARIUS**

# Your advantages with BRUMABA

## By Professionals, for Professionals

Their experience and procedures from everyday life in the operating theater are an important part of our product development. We have been working closely with surgeons and OR staff since 1980.

This is what makes BRUMABA so valuable for you.

- The products are specially designed to meet the requirements in Aesthetic, Plastic, Reconstructive Surgery and Dermatology.
- You and your team work with our operating tables and chairs in an ergonomic position at all times.
- Your patients are positioned professionally and comfortably during the entire surgical procedure - regardless of their size and weight.



Operating table  
**VARIUS**

**Flexible to the  
finest detail**





# The art of surgery

## The operating theater is the surgeon's stage

Every movement counts in the operating theater: Precision and concentration are crucial for the success of the operation and the well-being of the patient. Your surgical team works hand in hand, supported by essential instruments such as your BRUMABA operating table and operating chair. These help you to deliver your best performance at all times - and in the long term.

Thanks to the slim BRUMABA design, you enjoy optimum access to all areas of the patient while operating. This allows you to work closer to the operating field and adopt ergonomically ideal positions. Even areas, that are hard to reach, can be accessed effortlessly thanks to the wide range of adjustment options. This allows you to operate precisely in a comfortable position.

BRUMABA enables you to work ergonomically and perform precise operations. Experience maximum safety and efficiency for your everyday procedures with our solutions.



Ideal position for  
facelift



**Maximum  
flexibility:**  
One operating  
table for  
all procedures



Hand surgery in  
natural situation



Beach chair position  
for breast surgery



Lateral support



Lithotomy position



Prone position



# Operating table **VARIUS**

## All in One - Our masterpiece

### Standout features

- Two separately adjustable and spreadable leg section, for optimum access between the legs
- Airplane function: Tilt the patient towards you to work without straining your back
- Memory function: save up to 8 frequently used positions for efficient surgical procedures
- Premium, stable and durable stainless steel construction for highest quality "Made in Germany"
- Extra thick padding for even better comfort during long procedures
- 1 - 2 weeks of battery operation per charge guarantees reliable OR operation
- Three-dimensionally adjustable head rest
- Extensive height adjustment with DT double telescopic column (625 - 1055 mm)
- Up to 300 kg patient weight with full functionality



Slim design  
and legroom

### Optional equipment

- Our recommendation: 4-EASY-DRIVE for easy adaptation to the room situation and cleaning
- Foot switch for sterile, intraoperative adjustment



Foot switch  
for convenient  
adjustment



Scan the QR code  
and experience VARIUS live  
in the product video





Airplane function  
with 17°  
Rotation angle



Ergonomic work  
between the legs, e.g.  
liposuction, intimate  
surgery



## Technical data

Total length	1,770 - 1,940 mm
Total width	555 mm
Net weight VARIUS	205 kg
Net weight VARIUS N	195 kg
Max. table load	300 kg
Height (with padding)	625 - 1055 mm
Head section inclination	31° to -44°
Back section inclination	61° to -14°
Seat tilt	43° to -31°
Inclination foot section	30° to -90°
Axial tilt left/right	17°

**ART. NO. O2.02.31 VARIUS**  
**ART. NO. O2.12.31 VARIUS N**

# Operating table **PRIMUS**

## A unique operating experience

### Standout features

- The airplane function enables perfect ergonomic positioning with longitudinal axis tilting
- Our three-dimensionally adjustable headrest elevates your surgical precision and comfort
- Crafted from high-quality stainless steel, ensuring stability and durability "Made in Germany"
- Flexibly adjustable armrests for relaxed positioning
- Extra thick cushions maximize patient comfort during long procedures
- Memory function allowing you to save up to 8 frequently used lying positions, improving operational efficiency
- Long-lasting battery life of 1-2 weeks, guaranteeing reliable performance
- Expanded height adjustment with DT double telescopic column (540 - 970 mm)
- Up to 300 kg / 660 lbs patient weight with full performance

### Optional features

- Our recommendation: 2-GUIDE-DRIVE with two directional wheels for easy adjustment and cleaning in the operation room
- Foot switch feature for sterile and convenient intraoperative adjustments



Scan the QR code  
and experience PRIMUS  
live in the product video



Head shell for  
e.g. rhinoplasty



Wireless battery  
operation  
for up to two weeks





High-quality stainless steel construction



Optimum ergonomics



## Technical data

Total length	1,770 -1,940 mm
Overall width	555 mm
Deadweight	185 kg
Max. table load	300 kg
Height (incl. cushion)	540 - 970 mm
Inclination head section	31° to -44°
Inclination back section	72° to -17°
Inclination seat section	32° to -11°
Inclination foot section	30° to -90°

**ART. NO. O2.01.31 PRIMUS**  
**ART. NO. O2.11.31 PRIMUS N**

# OR table **GENIUS**

## Excellent down to the smallest detail

### Standout features

- Three-dimensionally adjustable head shell
- High-quality, stable and durable stainless steel construction for the highest quality "Made in Germany"
- Extra thick padding for even better patient comfort during long procedures
- Memory function: storage of up to 8 frequently used positions for efficient surgical procedures
- 1 - 2 weeks of battery operation per charge guarantees reliable OR operation
- Height adjustment with ST telescopic column (590 - 880 mm)
- Up to 250 kg patient weight with full functionality

### Optional equipment

- More extensive height adjustment with DT double telescopic column (530 - 960 mm)
- Foot switch for sterile, intraoperative adjustment
- Mounting rail on foot and head section for attaching accessories (ART. NO. Z1.151 Mounting rail 250 mm, Head/seat section)

ART. NO. G4.01.23 GENIUS DTA  
ART. NO. G4.11.23 GENIUS DTN  
ART. NO. G4.01.12 GENIUS STA  
ART. NO. G4.11.12 GENIUS STN



Panel and rail construction




Wireless battery operation  
for up to two weeks

Optional backrest  
upholstery  
with integrated armrests



The operating table shown contains special equipment.

Convenient storage



## Technical data

Total length	1,803.5 - 1,991.5 mm
Total width	565 mm
DTA Net weight	160 kg
Net weight STA	140 kg
Max. Load capacity of the table	250 kg
Height DTA (incl. cushion)	530 - 960 mm
STA height (incl. cushion)	590 - 880 mm
Head section inclination	30° to -45°
Back section inclination ST	80° to -22°
Back section inclination DT	80° to -11°
Inclination seat section ST	39° to -9°
Inclination seat section DT	35 to -9°
Inclination foot section	28° to -90°



# OR table **VENUS:**

## 2-in-1: treatment table and operating table

### Standout features

- Magnetic neck pad and integrated nose slot for perfect tummy and back support
- High-quality, stable and durable stainless steel construction for the highest quality "Made in Germany"
- Extra thick padding for even better patient comfort during long procedures
- Memory function: storage of up to 8 frequently used positions for efficient OR procedures
- 1 - 2 weeks of battery operation per charge guarantees reliable OR operation
- Height adjustment with ST telescopic column (630 - 910 mm)
- Up to 250 kg patient weight with full functionality



### Optional equipment:

- More extensive height adjustment with DT double telescopic column (530 - 960 mm)
- Foot switch for sterile, intraoperative adjustment
- Mounting rail set for attaching accessories (ART. NO. Z1.304)

ART. NO. V1.01.12 VENUS STA

ART. NO. V1.01.23 VENUS DTA

ART. NO. V1.11.12 VENUS STN

ART. NO. V1.11.23 VENUS DTN



Prone position



Ergonomic positioning for surgery and treatment



Mounting rail set VENUS  
(6 pcs.) for additional  
accessories

**ART.NO. Z1.304**

## Technical data

Total length	1,963 mm
Total width	660 mm
DTA Net weight	145 kg
Dead weight DTN	165 kg
Net weight STA	145 kg
Net weight STN	125 kg
Max. Load	250 kg
Height (with padding)	630-1,060 mm
Height DTA (incl. cushion)	630-910 mm
Upholstery height	80 mm

# Trolleys for operating tables

## For maximum mobility and stability



### Fixed base plate

Stable and Stationary: Ensuring the table remains securely in one place.

ART. NO. F1.101 PRIMUS

ART. NO. F1.103 GENIUS



### 2-GUIDE-DRIVE – GENIUS/PRIMUS

Easily move or lock the operating table with one lever. Features two directional and swivel castors, as well as a lift-brake system.

ART. NO. F1.203

### OUR RECOMMENDATION



### 4-EASY-DRIVE

- 4 swivel wheels with connecting axle
- Lift-brake system for extremely smooth standing
- Intelligent directional guidance
- Maximum legroom
- Lightweight 120 mm wheels

ART. NO. F1.403



### 5-WHEEL DRIVE

- Selective 5th wheel for easy, precise movement in long corridors and around tight corners
- Lift-brake system for maximum stability
- Maximum legroom
- 120 mm easy-running wheels

ART. NO. F1.504

# Controls and charging systems for operating tables



## Foot switch type IV with retaining bracket

- Type IV foot switch with 8 memory positions

**ART. NO. Z1.185**



## Foot control with holding bracket IV

For positioning the foot control on the chassis 4-EASY-DRIVE, 5-WHEEL-DRIVE.

**ART. NO. Z1. 274**



## Manual control - VARIUS/PRIMUS

- Included as standard
- With 8 memory positions

**ART. NO. Z1.118**



## Integrated charging system - PRIMUS/VARIUS

Charge your operating table once a week by plugging a cable directly into.

**ART. NO. X.000269**



## Interchangeable battery system with external charging station - PRIMUS

Replacement battery for smooth non-stop operation, particularly suitable in large operating rooms with changing OR staff to ensure additional operational reliability.

**ART. NO. Z1.134** Charging station  
**ART. NO. Z1.017** Battery 24 V  
(one battery already included)



## Rechargeable battery version 24 V – GENIUS

For independent operation of your operating table. Running time during normal operation of up to 2 weeks without recharging.

**ART. NO. G4.0X.XX**



## Power supply 110-230 V

For stationary use, the mains operation is suitable for direct connection to a socket.

**ART. NO. 02.11.31** PRIMUS N  
**ART. NO. G4.11.23** GENIUS DTN  
**ART. NO. G411.12** GENIUS STN

# Accessories for operating tables

## Head positioning



**Horseshoe-shaped  
headrest, small**

200x210x80 mm

**ART. NO. Z1.043A**



**Head plate, small**

320x240 mm

**ART. NO. Z1.109**

**Head plate, large**

400x300 mm

**ART. NO. Z1.003**



**Hair transplant  
support "neck support"**

**ART. NO. Z1.032**





**Armrest OP PRO**  
with fixing strap and  
length adjustment  
**ART. NO. Z1.331**

**Leg holder - OP,**  
with 2 ball-and-socket  
joints (1 pair)  
**ART. NO. Z1.343**



**Hand surgery table**  
with telescopic support  
**ART. NO. Z1.351**

**Mounting rail 250mm,**  
solid – **GENIUS/VENUS**  
Required for use of the  
leg holder or the hand  
surgery table with **GENIUS**  
and **VENUS**  
**ART.NO. X.000087**



## Accessories

For patient positioning like on cloud nine



**Back section -  
extension, narrow**

Included with VARIUS, PRIMUS  
and GENIUS



**Foot section extension /  
side extension, narrow - OP**

**ART. NO. Z1.267**  
PRIMUS/VARIUS

**ART. NO. Z1.263**  
GENIUS



**Arm support lateral,  
with fastening belt**

**ART. NO. Z1.002**



**Standing plate -  
VARIUS/VENUS (1 pair)**

**ART. NO. Z1.057**

**Lateral positioning  
support, deep**

180 x 300 x 50 mm

**ART. NO. Z1.110**



**Standing plate -  
VARIUS/VENUS (1 pair)**  
**ART. NO. Z1.112**




**Mattress viscoelastic,  
short with  
gynecological cutout**  
ART. NO. Z1.213



**Mattress  
viscoelastic**  
ART. NO. Z1.171



**Half roll cushion**  
ART. NO. Z1.100 (small)  
ART. NO. Z1.101 (large)



**Abdominal positioning  
cushion „chest“**  
ART. NO. Z1.345

**Abdominal positioning  
cushion „hips“**  
ART. NO. Z1.346



**Lordosis support cushion**  
ART. NO. Z1.140



**Stomach positioning  
cushion, trapezoidal  
(1 piece)**  
ART.NO. Z1.098

# Instrument tables

## Instrument table on the column

Plate size 600 x 400 mm

ART. NO. Z1.275

Plate size 700 x 500 mm

ART. NO. Z1.278



Infusion stand on the column  
ART. NO. Z1.280



Aluminum shelf for  
Instrument table  
ART. NO. Z1.245

Mobile  
Instrument table  
Plate size 600x400 mm  
ART. NO. Z1.172  
  
Plate size 700x500 mm  
ART. NO. Z1.282



## Remote control holders



Remote control holder  
on  
infusion stand  
ART. NO. Z1.243



Remote control holder  
FLEX, hinged  
ART. NO. Z1.313



Remote control  
holder, short  
ART. NO. Z1.307



Magnetic plate for  
remote control  
ART. NO. Z1.045





# Anesthesia accessories



**Anesthetic tube holder  
FLEX, foldable  
ART. NO. Z1.334**

**Infusion stand on the  
column, hinged  
ART. NO. Z1.277**



**Infusion stand with  
ball-and-socket joint  
ART. NO. Z1.038**



**Anesthetic frame,  
universal  
ART. NO. Z1.015**

# Protection, hygiene & assembly



**Mounting bracket**

ART. NO. Z1.004



**EasyCover**  
(150 pieces) disposable  
table cover

ART. NO. Z1.138



**Paper roll holder**

ART. NO. Z1.294 VARIUS

ART. NO. Z1.324 VENUS

ART. NO. Z1.051 PRIMUS/GENIUS



**Mounting glider with clamp**

ART. NO. Z1.093

**Protective Cover  
for foot section**

ART. NO. Z1.332

ART. NO. Z1.148A  
VARIUS



**Mounting adapter,  
25x10 mm**

For attaching  
accessories from other  
manufacturers

ART. NO. Z1.047



**Mounting rail  
– GENIUS/VENUS**

Included as standard  
with GENIUS

250 mm Head and seat section  
**ART.NR. Z1.151**

340 mm back and foot section  
**ART.NR. Z1.079**

# Patient recovery chair **CADDY**

## Design and comfort in a wake-up table

The CADDY is a 4-in-1 solution that functions as a chair, couch, transporter, and recovery chair. It is the perfect all-rounder that guarantees safe and comfortable patient positioning, coupled with easy transportation.

### Standout features

- Robust frame crafted from high-quality stainless steel
- Two multifunctional handles on the upper edge of the backrest for ergonomic pushing and easy adjustment of the reclining surfaces
- Comfortable cushions
- Physiologically advantageous lying positions ensure stable circulation for your patients
- Easy transport thanks to particularly smooth-running wheels with steering brakes
- Fold-up armrests ensure convenient entry and better access to the patient
- Choice of light or dark wood

### Optional accessories

- Armrest cushion with ball joint (ART. NR. Z1.070)
- Bodybelt, silicone (ART. NR. Z1.078)
- Caddy - Side protection (1 pair)
- Infusion stand (ART. NR. Z1.068)
- CADDY - wooden optic "beech light" standard (ART. NO. C1.XX1 )
- CADDY - wooden optic "beech dark" (ART. NO. C1.XX1, surcharge)



Scan the QR code and experience CADDY live in the product video



Side protection

Armrests automatically adjusted to reclining position (foldable)

Lying surface infinitely adjustable



## Technical data

Weight	67 kg
Max. dimensions (W x L x H) „Chair“	660 x 960 x 1440 mm
Max. dimensions (W x L x H) „Lounge“	660 x 1900 x 560 mm
Seat dimensions (W x L)	460 x 490 mm
Seat height	560 mm
Max. load	180 kg

ART. NO. C1.02.13

## Comfortable in any position



Chair position



Chaiselongue position



Shock position

# BRUMABA OR stools and chairs

## Invest in your health

BRUMABA operating stools and chairs can be flexibly adapted to your individual body dimensions and are the ideal complement to our operating tables.

Back protection minimizes physical strain during surgery and assisting at – a sustainable investment in your own health and that of your team.

### The Standout features of BRUMABA operating chairs

- Optimum back support thanks to infinitely adjustable back section
- Precisely configurable armrests for perfect support during the operation
- Patented ergonomic seat surfaces for a consistently secure and comfortable fit
- Individually adjustable tilt angle of the seat for healthy bloodcirculation in your legs
- Adjustable seat height for an angle of  $>90^\circ$  between your upper body and legs







Stool model  
**COMFORT**



Backrest adjustment



Height adjustment with  
foot release



Optimal  
Back support

Ergonomic  
Seat surface

Vein-friendly  
Sitting

Operating chair  
**BALANCE**

# Microsurgery chair **BALANCE ADVANCE**

## The comfortable and ergonomic solution for precision work

In close cooperation with microsurgeons, we have developed the BALANCE ADVANCE operating chair specifically for the requirements of microsurgery. The result: comfortable surgery with a perfect working position.





## Experience unmatched comfort and ergonomics

- Three-directional adjustable armrests with lockable feature support shoulder and arm muscles
- Comfortable work position even during extended procedures
- Ergonomically shaped seat and back with patented seat inclination and high-quality cushions ensure a healthy sitting posture
- Easy height adjustment with foot release (gas spring)
- Comfortable chassis with antistatic wheels
- High-frequency seamless cushions and closed wheels for hygienic cleaning

**BALANCE ADVANCE F**  
**ART. NO. S3.10.13.5310**

**BALANCE ADVANCE F&B (with brake)**  
**ART. NO. S3.11.13.5310**

**BALANCE ADVANCE SLIM F**  
**ART. NO. S3.10.23.641.0**

**BALANCE ADVANCE SLIM F&B (with brake)**  
**ART. NO. S3.12.23.641.0**

## Technical data

Min. height	555 mm
Max. height	745 mm
Chair cross diameter	610 mm
Seat cushion length	460 mm
Seat cushion width	420 mm

### BALANCE ADVANCE SLIM

Min. height	490 mm
Max. height	630 mm
Seat cushion length	400 mm
Seat cushion width	360 mm

# Operating chair **BALANCE**

## Achieving perfect balance: Form, function, and ergonomics

### **BALANCE F**

#### **Standout features**

- Body-adaptive seat and backrest for ergonomic support
- High-quality cushions for comfortable and healthy sitting
- The ergonomic seat design promotes healthy blood circulation in your legs, enabling you to sit comfortably for extended periods of time
- The cushions are seamlessly designed and the wheels are enclosed, allowing effortless cleaning and ensuring impeccable hygiene
- Optional column elevation (100 mm) (ART. NO. E.001558)



Optional with brake



Height adjustment  
with foot release,  
ideal for sterile work

#### **Technical data**

Min. height	500 mm
Max. height	695 mm
Chair base diameter	610 mm
Seat padding length	460 mm
Seat padding width	420 mm

**BALANCE F**  
**ART. NO. S3.10.13.501.0**

**BALANCE F&B (with brake)**  
**ART. NO. S3.11.13.501.0**



## **BALANCE PLUS F**

The BALANCE PLUS features adjustable armrests for optimal comfort and hygiene.

- Height and length of armrests can be customized to prevent tired arms
- Made of high-quality manmade material, the armrests are easy to clean and maintain high standards of hygiene
- Plain black design complements any color version of the chair

**BALANCE PLUS F**  
**ART. NO. S3.10.13.511.0**

**BALANCE PLUS F&B** (with brake)  
**ART. NO. S3.11.13.511.0**



## **BALANCE SUPREME F**

The BALANCE SUPREME is equipped with premium armrests that are padded and multifunctional.

- The armrests are adjustable in height, length, and can be moved to the side, providing a comfortable and customizable sitting experience
- Select the color of the armrest cushions to match the seat surface

**BALANCE SUPREME F**  
**ART. NO. S3.10.13.521.0**

**BALANCE SUPREME F & B** (with brake)  
**ART. NO. S3.11.13.521.0**



# Operating stool **COMFORT**

## The allrounder

Designed to enable medical staff to work with speed and agility, this maneuverable stool offers maximum flexibility for seamless use. With its practical foot-operated height adjustment, users can easily achieve the desired height while keeping their hands free.

### Standout features

- Seamless seat cushions for enhanced comfort and easy cleaning, ensuring reliable hygiene
- Large, enclosed, antistatic wheels that facilitate fast and agile working, while maintaining optimal hygiene
- Optional locking brake for a secure and stable working position
- Optional 100 mm column elevation for increased versatility and customizable height adjustment (ART. NO. E.001558)



### Technical data

Min. height	540 mm
Max. height	735 mm
Chair cross diameter	610 mm
Seat diameter	400 mm
Flexible height adjustment	510 - 710 mm



Height adjustment with foot release, ideal for sterile work



Available with COMFORT SWING seat cushion



## COMFORT ASSIST H

With 360° swivelling backrest for ergonomic support in a variety of sitting positions.

Only available with height adjustment via manual release.

**ART. NO. H3.10.41.103.0**



Simple height adjustment of the backrest

## Available COMFORT models

### COMFORT F

Height adjustment by foot release

**ART. NO. H3.10.11.100.0**

### COMFORT F/B

Height adjustment by foot release, with brake

**ART. NO. H3.11.11.100.0**



Optionally with brake

### COMFORT SWING F

With COMFORT SWING seat cushion, height adjustment by foot release

**ART. NO. H3.10.12.100.0**

### COMFORT SWING F/B

With COMFORT SWING seat cushion, height adjustment by foot release, with brake

**ART. NO. H3.11.12.100.0**



Height adjustment via hand release

### COMFORT H

Height adjustment by hand release

**ART. NO. H3.10.31.100.0**

### COMFORT SWING H

With COMFORT SWING seat cushion, height adjustment via hand release

**ART. NO. H3.10.32.100.0**

# BRUMABA products at a glance



**PRIMUS**



**VARIUS**



**GENIUS**



**VENUS**

# Multifunctional operating tables and operating chairs for Aesthetic, Plastic, Reconstructive Surgery and Dermatology



**CADDY**



**BALANCE  
PLUS**



**BALANCE  
SUPREME**



**BALANCE**



**BALANCE  
ADVANCE**



**COMFORT**



**COMFORT  
ASSIST**

# High-quality cushions

Maximize your comfort and hygiene

## Viscoelastic comfort cushions

Experience the ultimate in pressure relief and comfort with our viscoelastic cushions, designed to effectively prevent bedsores. From start to finish, we prioritize your comfort and that of your patients, ensuring a seamless and relaxing procedure.



## Elevate your hygiene standards to new heights

- The cushions can be easily removed for fast cleaning and disinfection
- Seamless and waterproof design to provide effective protection against bacteria
- High-quality durable artificial leather
- Our cushions are specially designed to meet all the hygiene requirements for operating rooms, ensuring the utmost safety and cleanliness for patients and staff
- Available for individual reordering, providing the ultimate flexibility for your seating needs



# Unleash your creativity

## 35 vibrant cushion colors

Our cushion colors are available for both operating tables and stools/chairs, allowing for seamless integration of your BRUMABA products into your operating room environment.



Scan QR code  
for further  
information



Color swatches may differ from the original color tone.



**BRUMABA GmbH & Co. KG**

Bürgermeister-Graf-Ring 17  
D-82538 Geretsried

Phone: +49 (0) 8171 / 2672 - 0

Fax: +49 (0) 8171 / 2672 - 10

Email: [info@brumaba.de](mailto:info@brumaba.de)

Subject to change without notice.

© BRUMABA GMBH & CO. KG

KA.0015-2024-10-22-EN

**[www.brumaba.com](http://www.brumaba.com)**

**MADE IN GERMANY**



**SINCE 1980**

**BRUMABA**  
OPERATING TABLE SYSTEMS