

DESIGNED, ENGINEERED AND MANUFACTURED  
SINCE 1980



**MADE IN GERMANY**



# VARIUS

A universal all-rounder  
among operating tables

SUITABLE  
SPECIALIST FIELDS:

**EYE**  
SERIES

OPHTHALMOLOGY

**OMS**  
SERIES

ORAL SURGERY  
MAXILLOFACIAL  
SURGERY

**OUTPATIENT**  
SERIES

OUTPATIENT SURGERY  
ORTHOPAEDICS  
UROLOGY  
GYNAECOLOGY  
PROCTOLOGY

**BEAUTY**  
SERIES

AESTHETIC SURGERY  
PLASTIC SURGERY  
HAIR TRANSPLANTATION

**SKIN**  
SERIES

DERMATOLOGY



OPERATING TABLE  
ART. NR. 02.02.31

# VARIUS

## A universal all-rounder among operating tables

The name says it all.  
With its high functionality, VARIUS covers almost all areas of specialty!

### EQUIPMENT FEATURES

- » Dual telescopic column 610-1040 mm
- » Longitudinal axis tilting
- » Memory function with 8 programmes
- » Split leg plates
- » Telescopic head section with head calvaria
- » Back section extension for length adjustment
- » Lever and steering castor undercarriage
- » 24V DC battery pack operation
- » Hand-held remote control
- » Stainless steel construction

### UNDERCARRIAGE (OPTIONAL)



#### 2-GUIDE-DRIVE "OP" - VARIUS

- » 2 Steering wheels with lever mechanism
- » 2 Directional rolling wheels, with connecting axle  
(L: 995mm - W: 508mm - Wheels: 120mm & 80mm)



#### 4-EASY-DRIVE - VARIUS

- » 4 Directional rolling castors, with connecting axle and lift-brake-system
- » Intelligent steering mechanism  
(L: 1065mm - W: 600mm - Wheels: 120mm)




PRONE POSITION



HEAD SURGERY



**BEACHCHAIR POSITION**



**KNEE ARTHROSCOPY**



**LATERAL POSITION**



**SHOULDER ARTHROSCOPY**



**DORSOSACRAL POSITION**



**X-RAY POSITION**




**KNEE-ELLBOW POSITION**



**SURGERY WHILE STANDING**


# Colour selection




## ELECTRIC CONDUCTIVITY

This pad meets all necessary standards for the use in an operating suite.


NO. 56001




F.4340-10120 / 56029



F.4340-20130 / 56079



F.4340-07422 / 56009




F.4340-07432 / 56019



F.4340-20150



F.4340-20162



F.4340-20161




F.4340-20156



F.4340-07496 / 56068



F.4340-20154



F.4340-07445 / 56040



F.4340-00002 / 56024




F.4340-5347



F.4340-50454



F.4340-20157



F.4340-07459 / 56037



F.4340-10295 / 56066



F.4340-07436 / 56023



F.4340-07480 / 56056



F.4340-50621



F.4340-20149



F.4340-07478 / 56053



F.4340-07479 / 56054



F.4340-50618



F.4340-07485 / 56048



F.4340-50620




F.4340-07476 / 56051



F.4340-07477 / 56052



F.4340-20127 / 56080



F.4340-50619



F.4340-50622



F.4340-50623



F.4340-5015



F.4340-20129 / 56083



**NEW:** Also available in VISCOELASTIC form.

DESIGN, CONSTRUCTION  
& MANUFACTURING

# MADE IN GERMANY

SINCE 1980

## BRUMABA GMBH & CO. KG

Bgm.-Finsterwalder-Ring 21  
D-82515 Wolfratshausen  
Germany

Telephone: +49 (0) 8171 / 2672 - 0  
Telefax: +49 (0) 8171 / 2672 - 10

Email: [info@brumaba.de](mailto:info@brumaba.de)  
Internet: [www.brumaba.de](http://www.brumaba.de)

All changes reserved | Edition 2016  
© BRUMABA GMBH & CO. KG

Printnr.: 2016\_02 000014

On request we will send you our binding colour and material samples.